

PRIMARY ANTIBODY CARTRIDGE (PAC) PROGRAMMER

REF 014-2524-000

INSTRUCTIONS FOR USE
PAGE 1 of 1

INTENDED USE

IVD For *in vitro* diagnostic use.

The PAC Programmer is intended for laboratory use to program the Smart Tags of a User Fillable Primary Antibody Cartridge (PAC). The PAC Programmer is intended for use exclusively with the Celerus Wave[®] RPD System.

INSTRUCTIONS FOR USE

Refer to the Celerus Wave[®] RPD Operator's Manual for additional information including Materials Required but Not Provided; Storage; Staining Procedure; Troubleshooting; Interpretation of Staining; and General Limitations.

MATERIALS REQUIRED BUT NOT PROVIDED

Refer to the Celerus Wave[®] RPD Operator's Manual for additional information about Materials Required but Not Provided.

STORAGE AND HANDLING

Refer to the Celerus Wave[®] RPD Operator's Manual for additional information about Storage and Handling.

WASTE DISPOSAL











The PAC Programmer should be returned to Celerus Diagnostics for environmentally safe disposal.

PRECAUTIONS

For professional users.

Refer to the Celerus Wave[®] RPD Operator's Manual for additional precautions.

EXPLANATION OF SYMBOLS

| | | | | | |
|---|---|---|--|---|---|
| REF Catalog Number |  Temperature Limitations | LOT Batch Code |  Sufficient for <n> Tests |  Harmful, Irritant |  Warning |
| IVD In Vitro Diagnostic Medical Device |  Use By |  European Conformity | EC REP EC Representative |  Manufacturer |  Refer to Instructions for Use |

