

Technical Specifications

Dimensions	110 x 60 x 70 cm 43 x 23.5 x 27.5 in (W x H x D)
Weight	57 kg; 126 lbs
Electrical Requirements	100-240 VAC; 50-60 Hz; 10 Amps; 510 W
Fuse Requirements	None
Environmental	Operation: 5 °C-40 °C (41 °F-104 °F) at 5-80% relative humidity Storage: -20 °C-65 °C at 5-95% relative humidity
Total Slide Capacity	2-16 glass slides (US & International sizes)
Reagent Capacity	Up to 8 each 100-Test Antibody PAC and 100-Test LRM Detection Chemistries

Components included with the Celerus Wave™ System

Label Printer	Zebra Label Printer 3844-Z 19.8 x 17.3 x 23.8 cm 7.8 x 6.8 x 9.4 in (W x H x D)
Slide Labels	Heat and chemical resistant
Computer	Intel Celeron 2.66 GHz laptop
Printer	USB compatible letter size printer; All required cables included
Warranty	12 months parts and labor